

Product Name :
Suction Unit - HIVAC JR.**Product Code :**
NMSUC10HV5**Description :**

Suction Unit - HIVAC JR. Housing: MS powder coated on four castors Capacity: -650 mm Hg \pm 10 , 42 Ltrs/ min
Pump Type: Piston Pump Jars: 2 x 2 Ltrs Polycarbonate Jar Filter: Trap Filter Tubing: 10 mm ID x 2 mtr. (PVC)
Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 130Watt (110V on request)
Noise Level: 50 dB A \pm 3 Dimension & Weight: 39 x 31 x 62 cms, 15 Kg. Overflow Safety: Dual Mechanical

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers,](http://MedicalHospitalEquipmentsManufacturers.com)