

**Product Name :**  
Manual Suction Unit (Foot / Pedal Operated)

**Product Code :**  
NMSUC30HV3



**Description :**

Manual Suction Unit (Foot / Pedal Operated) Capacity: Foot Operated Vacuum pump giving vacuum more than 550mm ( + 50mm). Jars: Two Autoclavable PC jars with cap of 1 litre capacity. A suction cannula is also provided with the unit. The unit is fitted on a powder coated base

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)