

Product Name :
Suction Unit HI-VAC Plus (90 Ltrs./Min)**Product Code :**
NMSUM10HV1**Description :**

Suction Unit HI-VAC Plus Tender Suppliers Africa and Suction Unit HI-VAC Plus Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Suction Unit HI-VAC Plus Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Suction Unit HI-VAC Plus China, Suction Unit HI-VAC Plus Kenya.

Suction Unit HI-VAC Plus Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Suction Unit HI-VAC Plus Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the Grenadines

(Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point),
Suriname (Paramaribo)

Technical Specification :

Suction Unit HI-VAC Plus (90 Ltrs./Min) Housing: Plastic moulded on two wheels Capacity: -700 mm Hg \pm 10,
90 ltrs/ min Pump Type: Double Piston Pump Jars: 2 x 2.5 Ltrs wide mouthed (Polycarbonate) with changeover
lever. Provided with handle for hygiene while emptying/cleaning. Filter: Bacterial Filter Autoclavable/Reusable
Tubing: 10 mm ID x 2 mtr. (PVC) Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 /
60 Hz, 580Watt (110V on request) Noise Level: 60 dB A \pm 3 Dimension & Weight: 46 x 34 x 86 cms, 23 Kg.
Overflow Safety: Mechanical / Electronic Acoustic Alarm Accessories: Foot Switch

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)