

**Product Name :**  
Suction Unit MAGNAVAC**Product Code :**  
NMSUM10HV3**Description :**

Suction Unit MAGNAVAC (70 Ltr/Min) High Vacuum Suction Unit for LIPO SURGERY Housing: Stainless Steel on four castors

**Technical Specification :**

- High Vacuum Suction Unit for LIPO SURGERY
- **Housing:** Stainless Steel on four castors
- **Capacity:** -700 mm Hg  $\pm$  10 , 70 Ltrs/ min
- **Pump Type:** Double Rotary Vane Pump
- **Jars:** 2 x 2 Ltrs Polycarbonate Jar
- **Filter:** Trap Filter with inbuilt overflow protection.
- **Tubing:** 10 mm ID x 2 mtr. (PVC)
- **Vacuum Gauge:** 10 cm dia 0-760 mm Hg calibration.
- **Power:** 220 V AC, 50 / 60 Hz, 380Watt (110V on request)
- **Noise Level:** 50 dB A  $\pm$  3
- **Dimension & Weight:** 49 x 36 x 70 cms, 36 Kg.
- **Overflow Safety:** Dual Mechanical

---

Website: [SITE](#)

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**