

**Product Name :**  
Suction Unit - Stand Model**Product Code :**  
NMSUM10HV4**Technical Specification :**

Suction Unit - Stand Model Housing: ABS Plastic Moulded on Trolley with four castors. Capacity: -575 mm Hg  
Flow: 10, 45 Ltrs/ min Pump Type: Diaphragm Pump Jars: 2 x 2 Ltrs Polycarbonate Jar Filter: Bacterial Filter  
Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (Silicon) Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration.  
Power: 220 V AC, 50 / 60 Hz, 90Watt (110V on request) Noise Level: 55 dB A  
Dimension & Weight: 48 x 48 x 114 cms, 8 Kg + 9 Kg Trolley Overflow Safety: Mechanical / Electronic Acoustic Alarm Accessories: Foot Switch

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)