

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Hot Water Bottle Cover - Knitted, 100% Cotton

Product Code :
MRG20100030



Description :

Hot Water Bottle Cover - Knitted, 100% Cotton Knitted Hot Water Bottle Cover, 100% Cotton Available in any color Washing Instructions: - Machine wash - Do not bleach - Tumble dry - Cool iron if required

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE,](#)