

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Tourniquet - Latex

Product Code :
MRG20100057



Description :

Tourniquet - Latex Features: Tourniquet Latex Elastic Strip Tubular Reusable, Non-sterile. Length: 75-100cm Approx. Instructions for use: A tourniquet is used to facilitate the puncture of veins. Safety Process: This item must be cleaned and disinfected as often as necessary

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,