

Product Name :
Infant Resuscitation Kit (Infant)**Product Code :**
NMCBS26044**Description :**

Infant Resuscitation Kit (Infant) An Essential Carrying Kit for providing Effective Resuscitation to an Infant. (Available in a Carry Bag). Contents of the Kit Infant Resuscitator, Foot Suction, Laryngoscope, Airways. Infant Resuscitator: Rugged; 100% Autoclavable; Reusable. 250 ml silicone bellows. Non re breathing valve with 40 cm of H₂O pressure release. Mask sizes : No. 00 & 01. 360° swiveling patient connector with standard 15 mm / 22 mm diameter. Corrugated PVC oxygen reservoir. 1.5 mtr. PVC oxygen tubing. Foot Suction: Compact, light weight, easy to handle and operate. Durable rubber bellows. Long lasting stainless spring to provide friction free pumping. Completely autoclavable polycarbonate vacuum jar with lid (500 ml capacity). Scratch proof powder coated frame. Laryngoscope : Stainless steel straight blades of size 0 & 1 with penlight handle. Airways : Silicone, autoclavable & reusable in sizes 000, 00 & 0

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers,](http://MedicalHospitalEquipmentsManufacturers.com)