

Product Name :
Resuscitation Unit with Fixed Baby Cradle & Drawers with
Microprocessor based temperature controller with 3 modes
Skin/Air/Manual

Product Code :
NMCBS26046



Description :

Resuscitation Unit with Fixed Baby Cradle & Drawers with Microprocessor based temperature controller with 3 modes Skin/Air/Manual Technical Specification Heater Box Heating Element : Ceramic Heater 650W with parabolic reflector, protected by metallic grid The source box can be swivelled to take X-ray Halogen Examination Lamp (12 Volts, 50 Watts) Provides accurate assessment of Infant colour and sufficient illumination during procedures Baby Tray with mattress Transparent acrylic, drop down and lockable side panels Baby Tray with mattress, along with head up/down facility and X-ray cassette holder Mattress density approx 21-25Kg/m³ and with removable, washable, waterproof cover Acrylic X-ray tray with X-ray guide for positioning the cassette The design aspects also allow ease of dismantling baby care area during cleaning or fumigation Essential Attachments Lower Shelf and plastic moulded 3 Storage Drawer Side tray for Monitor Pole for Syringe Pump / Infusion Pump S.S. IV Rod 4 swivel castors (100mm), 2 castors with brake Oxygenation : Integrated support for 10L oxygen Cylinder Oxygen regulator with pressure gauge, flow meter and humidifier bottle along with silicone tubing for delivery of oxygen from humidifier bottle to patient Suction : Integrated Slow Suction unit with Vacuum Regulator, Vacuum Gauge and Autoclavable Fluid Jar 500ml (2 nos.) and Silicone tubing Resuscitation : Silicone manual Resuscitator size Infant (250ml) with 00 and 01 size face masks. Laryngoscope - Penlight handle with size 0 and 1 straight blade Silicone Airways size 000, 00 and

Website: Medical Hospital Equipments Manufacturers,

NAUGRA[®]

Medical Hospital Equipments Manufacturers