

**Product Name :**  
Resuscitation Unit with Fixed Baby Cradle & Drawers with  
Microprocessor based temperature controller with 2 modes  
Skin/Manual mode

**Product Code :**  
NMCBS26047



**Description :**

Resuscitation Unit with Fixed Baby Cradle & Drawers with Microprocessor based temperature controller with 2 modes Skin/Manual mode  
Technical Specification  
Heater Box Heating Element : Ceramic Heater 650W with parabolic reflector, protected by metallic grid  
The source box can be swivelled to take X-ray Halogen Examination Lamp (12 Volts, 50 Watts)  
Provides accurate assessment of Infant colour and sufficient illumination during procedures  
Baby Tray with mattress  
Transparent acrylic, drop down and lockable side panels  
Baby Tray with mattress, along with head up/down facility and X-ray cassette holder  
Mattress density approx 21-25Kg/m<sup>3</sup> and with removable, washable, waterproof cover  
Acrylic X-ray tray with X-ray guide for positioning the cassette  
The design aspects also allow ease of dismantling baby care area during cleaning or fumigation  
Essential Attachments  
Lower Shelf and plastic moulded 3 Storage Drawer  
Side tray for Monitor Pole for Syringe Pump / Infusion Pump  
S.S. IV Rod 4 swivel castors (100mm), 2 castors with brake  
Oxygenation : Integrated support for 10L oxygen Cylinder  
Oxygen regulator with pressure gauge, flow meter and humidifier bottle along with silicone tubing for delivery of oxygen from humidifier bottle to patient  
Suction : Integrated Slow Suction unit with Vacuum Regulator, Vacuum Gauge and Autoclavable Fluid Jar 500ml (2 nos.) and Silicone tubing  
Resuscitation : Silicone manual Resuscitator size Infant (250ml) with 00 and 01 size face masks. Laryngoscope - Penlight handle with size 0 and 1 straight blade  
Silicone Airways size 000, 00 and

---

Website: Medical Hospital Equipments Manufacturers,

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**