

Product Name :
Infant Incubator Single Wall (Single Canopy) without drawers with Microprocessor based temperature controller with Skin / Air / Manual mode

Product Code :
NMCBS26052



Description :

Infant Incubator Single Wall (Single Canopy) without drawers with Microprocessor based temperature controller with Skin / Air / Manual mode
Technical Specification: Acrylic Single Wall Canopy with six port holes. Front Loading. Four Inlets for administration of O2, IV lines or Ventilator tubing. Sliding Baby Tray with External jerk free gear based Head Up / Down facility. X-Ray cassette sliding facility. Stainless Steel rust free Heater, blower sub assembly with changeable micro Air filters. Provision of an open shelf in the trolley. Humidity heater and Hygrometer

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,