

**Product Name :**  
Infant Radiant Warmer (ECO-BASIC)**Product Code :**  
NMCBS26058**Description :**

Infant Radiant Warmer (ECO-BASIC) Technical Specification Microprocessor based servo controlled unit allow two operation mode, automatic and manual mode. Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. Transparent Side panel. Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. A second thermal cut-out function for more safety. Skin temperature sensor failure function to avoid over temperature

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)