

Product Name :
Infant Radiant Warmer (ECO-TOP)**Product Code :**
NMCBS26059**Description :**

Infant Radiant Warmer (ECO-TOP) Technical Specification Microprocessor controlled unit allow three operational mode: Pre-warm, Manual or Automatic(Servo Controlled temperature System). Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. Transparent Side panel. Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. Heater head can be adjusted $\pm 90^\circ$ horizontally. Convenient X-ray cassette board under infant bed. Inclination of infant bed is adjustable. A second thermal cut-out function for more safety. Skin temperature sensor failure function to avoid over temperature

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,