

**Product Name :**  
Infant Radiant Warmer (STANDARD)**Product Code :**  
NMCBS26060**Description :**

Infant Radiant Warmer (STANDARD) Technical Specification Microprocessor controlled unit allow three operational mode: Prewarm, Manual or Automatic(Servo Controlled temperature System). Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. APGAR timer to record therapy time Convenient X-ray cassette under infant bed Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. Heater head can be adjusted  $\pm 90^\circ$  horizontally Inclination of infant bed is adjustable. Observe light angle is adjustable. Foldable transparent side panels  
Optional Parts 1. RS232 Connector 2. Drawer 3. Height Adjustable