

**Product Name :**  
Resuscitation Unit (PREMIUM)**Product Code :**  
NMCBS26064**Description :**

Resuscitation Unit (PREMIUM) Microprocessor controlled unit allow three operational mode, Pre-warm, Manual or Automatic (Servo Controlled temperature System). Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. APGAR timer to record therapy time Convenient X-ray cassette under infant bed Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. Heater head can be adjusted  $\pm 90^\circ$  horizontally Inclination of infant bed is adjustable. Observe light angle is adjustable Foldable transparent side panels. A second thermal cut-out function for more safety Skin temperature sensor failure function to avoid over temperature. Low pressure suction unit, oxygen supply system and LED phototherapy unit Optional Part: 1. RS232 connector 2. Drawer

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,