

**Product Name :**  
Macintosh Type Curved Laryngoscope Blades - Stainless Steel, POLISHED FINISH

**Product Code :**  
NMLYG10B07



**Description :**

Macintosh Type Curved Laryngoscope Blades - Stainless Steel, POLISHED FINISH Laryngoscope Blade SS. Blade Size 4. Polished. Autoclavable

**Technical Specification :**

Macintosh Laryngoscope Blades combine durable stainless steel with large diameter, high performance full integrated fiber optic bundles.

Usage: Hospital

Material: Stainless Steel

Automation Grade: Automatic

Product Type: Laryngoscope

Shape: Curved