

Product Name : Macintosh Type Curved Laryngoscope Blades - Stainless Steel, POLISHED FINISH	Product Code : NMLYG10B08
-------------------------------------------------------------------------------------------------------	-------------------------------------



Description :
Macintosh Type Curved Laryngoscope Blades - Stainless Steel, POLISHED FINISH Laryngoscope Blade SS. Blade Size 5. Polished. Autoclavable
Technical Specification :
Macintosh Laryngoscope Blades combine durable stainless steel with large diameter, high performance full integrated fiber optic bundles. Usage: Hospital Material: Stainless Steel Automation Grade: Automatic Product Type: Laryngoscope Shape: Curved