



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Miller Type Straight Laryngoscope Blades - Stainless Steel
(Satin/Dull Finish)

Product Code :
NMLYG10A12



Description :

Miller Type Straight Laryngoscope Blades - Stainless Steel (Satin/Dull Finish) Laryngoscope Blade SS. Blade Size 00. Matt

Technical Specification :

Combines the straight Jackson blade design with the curved distal tip Miller design to allow greater visibility and a 'direct line' approach to the trachea
Low profile flange reduces risk of damage to teeth
Crafted in surgical stainless steel
Satin matte finishes which eliminates reflection and glare
Fit any Standard Conventional handles



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,