

Product Name :
Miller Type Straight Laryngoscope Blades - Stainless Steel
(Satin/Dull Finish)

Product Code :
NMLYG10A16



Description :

Miller Type Straight Laryngoscope Blades - Stainless Steel (Satin/Dull Finish) Laryngoscope Blade SS. Blade Size 3. Matt

Technical Specification :

Combines the straight Jackson blade design with the curved distal tip Miller design to allow greater visibility and a 'direct line' approach to the trachea
Low profile flange reduces risk of damage to teeth
Crafted in surgical stainless steel
Satin matte finishes which eliminates reflection and glare
Fit any Standard Conventional handles