

**Product Name :**  
Macintosh Type Reusable Fibre Optic Laryngoscope Blade

**Product Code :**  
NMLYG10B15



**Description :**

Macintosh Type Reusable Fibre Optic Laryngoscope Blade Laryngoscope Blade - With Detachable Fibre Optic Tube Autoclavable

**Technical Specification :**

High grade of fibre-optics provides optimal light transfer without heat, reducing risk of thermal damage and maximising view of the larynx.

Made from high quality stainless steel, our blades give clinicians a complete choice for all patient intubation requirements.