

Product Name :
Miller Type Reusable Fibre Optic Laryngoscope Blade

Product Code :
NMLYG10A18



Description :

Miller Type Reusable Fibre Optic Laryngoscope Blade Laryngoscope Blade - With Detachable Fibre Optic Tube Autoclavable

Technical Specification :

High grade of fibre-optics provides optimal light transfer without heat, reducing risk of thermal damage and maximising view of the larynx.

Made from high quality stainless steel, our blades give clinicians a complete choice for all patient intubation requirements.