

<b>Product Name :</b> Laryngoscope set with Standard Illumination - in Satin/Dull Finish	<b>Product Code :</b> NMLYG10S08
---	-------------------------------------



<b>Description :</b>
Laryngoscope set with Standard Illumination - in Satin/Dull Finish The set of four satin blades, includes: Curved macintosh type, size: 1, 2 & 3. Straight miller type, size: 1. Spare bulbs - 4 pcs., & Stainless steel satin handle uses two alkaline batteries (medium size, 1.5v battery)
<b>Technical Specification :</b>
Smooth barrel knurling provides perfect grip without risk of tearing s gloves.
Stainless steel points increase life.
Fits all conventional blades made to this standard
SAFE and EASY TO USE design : Brushed satin finish reduces glare and fits to all conventional handles.
The LED light bulb in the blade can easily be changed.
Blade made from highest quality stainless steel.
Cool LED light protects delicate oral tissues

---

## Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

**NAUGRA<sup>®</sup>**