

Product Name :
Large Vaccine Carrier, Long Range**Product Code :**
VCS20100004**Description :**

Large Vaccine Carrier, Long Range Vaccine Storage Capacity as per Std. PQS load : 2.46 Litres. Gross Storage Capacity : 2.50 Litres. Weight Fully Loaded : 6.36 Kgs. Weight Empty (with empty ice pack) : 2.98 Kgs. External surface material : HDPE* Internal lining material : HIPS* Insulation material : CFC free Polyurethane Insulation thickness : 32/35 mm External Dimensions(LxWxH) : 27.0x27.0x32.0 cms. Internal Dimensions(LxWxH) : 12.2x12.2x19 cms. Vaccine Storage Dimension : 11.38x11.38x19.0 cms. Lid Type and Fixings : Removable Number of Ice Packs required : 4. Number of Ice Packs supplied : 4. Ice-pack type compliant with PQS EOS/IP01.1 : 0.6 L. Volume of Ice Pack : 0.60 Litres. Min. Cold life Req. : 30 hrs. Conforming to WHO PQS Ref. No. : E04/VC01.1. Cold Life without openings at +43°C : 50 hrs. 12 mins. Warm Life at -20°C : 15 hrs. 57 mins. Cool Life at +43°C : 9 hrs. 47 mins