


**Product Name :**  
Small Cold Box Short range**Product Code :**  
VCS20200005**Description :**

Small Cold Box Short range Vaccine Storage Capacity as per Std. PQS load : 16.0 Litres. Gross Storage Capacity : 20.80 Litres. Weight Fully Loaded : 34.42 Kgs. Weight Empty (with empty ice pack) : 14.30 Kgs. External surface material : LLDPE\* Internal lining material : LLDPE\* Insulation material : CFC free Polyurethane Insulation thickness : 65 mm External Dimensions(LxWxH) : 65.0x65.0x34.7 cms. Internal Dimensions(LxWxH) : 47.8x48.1x22.5 cms. Vaccine Storage Dimension : 40.1x40.6x12.78 cms. Lid Type and Fixings : Hinged Number of Ice Packs required : 32. Number of Ice Packs supplied : 32. Ice-pack type compliant with PQS E0S/IP01.1 : 0.4 L. Volume of Ice Pack : 0.4 Litres. Min. Cold life Req. : 48 hrs. Conforming to WHO PQS Ref. No. : E04/CB01.1. Cold Life without openings at +43°C : 80 hrs. 19 mins. Warm Life at -20°C : 27 hrs. 44 mins. Cool Life at +43°C : 18 hrs. 44 mins