

Product Name :
Long Cold Box Long range**Product Code :**
VCS20200006**Description :**

Long Cold Box Long range Vaccine Storage Capacity as per Std. PQS load : 18 Litres. Gross Storage Capacity : 21.70 Litres. Weight Fully Loaded : 46.14 Kgs. Weight Empty (with empty ice pack) : 21.34 Kgs. External surface material : LLDPE* Internal lining material : LLDPE* Insulation material : CFC free Polyurethane Insulation thickness : 100 mm External Dimensions(LxWxH) : 77.0x61.0x51.0 cms. Internal Dimensions(LxWxH) : 53.0x36.1x30.8x cms. Vaccine Storage Dimension : 45.0x29.4x16.4 cms. Lid Type and Fixings : Hinged Number of Ice Packs required : 44. Number of Ice Packs supplied : 44. Ice-pack type compliant with PQS E0S/IP01.1 : 0.4 L. Volume of Ice Pack : 0.4 Litres. Min. Cold life Req. : 96 hrs. Conforming to WHO PQS Ref. No. : E04/CB01.1. Cold Life without openings at +43°C : 147 hrs. 14 mins. Warm Life at -20°C : 52 hrs. 03 mins. Cool Life at +43°C : 30 hrs. 39 mins