

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Hand Held Suction Unit

Product Code :
NMEGY10025



Description :

Hand Held Suction Unit Handheld portable suction unit designed to provide quick and effective suction for multiple emergency use. Vacuum(max) : 450mmHg . Disposable Container Volume : 200 ml. Operating environmental temperature : - 20 °C to + 50 °C. Storage Environmental Temperature: - 40 °C to + 60 °C. Dimensions &Weight : 17 x 5 x 18 cm. 192 gm

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,