

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Oxygen Concentrators

Product Code :
Surgical Instruments0051



Description :

Oxygen Concentrators It apply zeolite molecular sieve as shift, adopt advanced physical theory(PSA) to separate oxygen and nitrogen in air. At the same time filter the bad substance in air, so as to get pure oxygen for medical use Oxygen Concentration : 90%. Oxygen flux: 1~5L/min. Output pressure: 0.04~0.05Mpa. Overrunning pressure alarm Noise: 55dBA. Operation noise: 49dbA. Power supply: AC220V±20V, 50Hz±1 Hz. Input power: 600VA

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)