

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Sphygmomanometer Aneroid Round Mobile

**Product Code :**  
NMDIG10015



**Description :**

Sphygmomanometer Aneroid Round Shaped Mobile, ABS Round dial face adjustable for easy reading, stand on aluminum tube with 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing. Sphygmomanometer Aneroid Round Mobile, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)