

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Sphygmomanometer Aneroid Round Mobile

Product Code :
NMDIG10015



Description :

Sphygmomanometer Aneroid Round Shaped Mobile, ABS Round dial face adjustable for easy reading, stand on aluminum tube with 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing. Sphygmomanometer Aneroid Round Mobile, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,