

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
U SPLINT

**Product Code :**  
NMBSK1018



**Description :**

U SPLINT Size: Short arm, Long arm, Short leg, Long leg. Features: Malleable, semi-rigid, light metal splint provides support in favourable anatomical position. Provides immobilization to the injured part. Splint is X - ray translucent. Applications: Emergency splinting of dislocated / fractured upper or lower limbs

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)