



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Macintosh Laryngoscope Blades

Product Code :
NMLYG10B02



Description :

Macintosh Type Curved Laryngoscope Blade Dull Finish Size 0, 1, 2, 3, 4, 5 Macintosh Laryngoscope Blades, Macintosh Laryngoscope Blades Exporters, Macintosh Laryngoscope Blades Manufacturers

Technical Specification :

Macintosh Laryngoscope Blades combine durable stainless steel with large diameter, high performance full integrated fiber optic bundles.
Designed with sterilization in mind, integrated system minimizes fluid entrapment. Epoxy sealed channel is resistant to temperatures up to 392°F.
Sterilize in Steam, EO and Gas Plasma Systems
Green system compliant
Reinforced ball-lock system
Flawless pinhole free Ni/Cr plating
Select quality 304 Stainless Steel



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,