

NAUGRA[®]

Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
ULTRA VIOLET LAMP

Product Code :
PEEHC10005



Description :

ULTRA VIOLET LAMP Floor stand model for Ultraviolet radiation treatments. UV bulb of 300w is fitted in metal casing with auto-matic cut off timer of 15 minutes. Adjustable arm permits wide range of adjustments. Mounted on heavy duty base fitted with ball bearing casters. Comes with suitable power cord

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)