

---

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
PARAFFIN WAX BATH

**Product Code :**  
PEEHC10011



**Description :**

PARAFFIN WAX BATH Internal wax tank made out of stainless steel sheet. The top has a aluminum cover and the edges of the bath has a laminated wooden rim, works on electricity 220-240 V. Provided with special wax immersion type heater and thermostat to control temperature. Externally finished oven baked. Wax Chamber size 40cm x 20cm x 20cm., Comes with Three kg wax

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](http://SITE)