

**Product Name :**  
ARM EXERCISER

**Product Code :**  
PEEEE10035



**Description :**

ARM EXERCISER All chrome plated, tubular steel construction, mounted on wooden platform and fitted with painted wall bracket at the top for wall mounting. Comes with adjustable swivel pulley block, a set of Ten locking weights and two handles with nylon rope, suitable for all shoulder, elbow & chest exercises

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,