

Product Name :
Hydrotherapy Tank**Product Code :**
PEEHT10001**Description :**

Hydrotherapy Tank Full body T shape hydrotherapy tank for hydro-massage. The heavy gauge stainless steel tank, rests on 40cm (16") high tubular stand re-inforced with cross bars. Overall size of the tank is 190cm(75") x 185cm (73") x 93cm (37"). Width of the narrow portion of the tank is 91cm (36") and depth inside the tank is 48cm (19"). Comes with Two motorized electric turbine ejectors & aerators with turbine carriages and elevators, Digital Thermometer, adjustable Head-Rest, Body Stretcher and Sitting seat fitted with nylon cloth. Tank is fitted with Two water inlets, One overflow and one Drain fittings. Tank is finished bright inside & painted outside

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,