

Product Name :
WHIRLPOOL BATH**Product Code :**
PEEHT10003**Description :**

WHIRLPOOL BATH Heavy gauge stainless steel tank sturdily fabricated and bright finished. Provides effective hydro-massage of arms, hips, legs, knees and feet. Mounted on heavy duty 10cm dia four ball bearing castors. Comes complete with removable inside aluminium Seat & Arm rest. Unit is fitted with a motorized turbine ejector and aerator on a spring loaded turbine elevator, digital thermometer, thermostat & immersion heater. Provided with drain pipe for emptying. Inside measurements of tank 90cm long x 50cm wide x 70cm deep. Inside measurements of tank 90cm long x 50cm wide x 70cm deep

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,