



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Dental Chair Programmable (Conventional Underarm  
Delivery Unit)

**Product Code :**  
Dental0002



**Description :**

Dental Chair Programmable (Conventional Underarm Delivery Unit) Electrically operated Dental chair unit with Up / Down / Back/ Forward. Auto return to zero position. Feathers touch control on delivery unit & chair backrest. Unit Consisting Followings: Air-Rotor terminal – 1 No. Micro - Motor terminal control – 1 No. Three way syringe -1 No. Two intensity operating light fitted with imported reflector. Water connection for spittoon and tumbler through electrically operated solenoid valve. Arm rest both sides for better patient comfort. Articulated adjustable headrest with facility to convert into pediatric headrest. Stainless Steel instrument Tray with hanging Straight Tubing. Dr. Stool – Pneumatic. X-Ray viewer – LED based. Seamless Fully laminated imported upholstery in attractive colors. Multifunctional foot control. Single piece Ceramic Spittoon top. Wide instruments tray. 2 Actuators (Noise, Vibration & Jerrck Free). 3 Side Operated Chair



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,