

**Product Name :**  
Ultra Sonic Scaler (With 3 Tips)

**Product Code :**  
DSUS0002



**Description :**

Ultra Sonic Scaler (With 3 Tips) Special self- contained water supply system makes it portable. Enable to use dental solutions for different treatment. Multifunction for general scaling, Endodontic, Periodontal, and restorative treatments. Solution bottle can be change within 3 seconds. Compatible with EMS tips. Built- in water-cooling system extend PCB lifetime. Auto-boost function provides sufficient power to remove heavy calculus. Top design: all surfaces are easy for clean and maintenance. Autoclavable hand piece can be autoclaved up to 135½C. No Compressor & No additional water system required