



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Dental H.P Oil Free & Noise Less Air Compressor (Fully Imported)

Product Code :
DentalAir0002



Description :

Dental H.P Oil Free & Noise Less Air Compressor (Fully Imported) Heavy-duty air reservoir tank made from thick M.S. sheet with Air & Tank Pressure gauge. Automatic cut off. Switch, safety valve, ON/OFF valve, Drainage valve etc. Moisturizer Filter With Regulator. Tank capacity 50 Lts. V/Hz 200V/50Hz. Easy handling



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,