

Product Name :
Light Cure Unit With Cord Digital LCD Display (LED based)**Product Code :**
DLCU0001**Description :**

Light Cure Unit With Cord Digital LCD Display (LED based) Tender Suppliers Africa and Light Cure Unit With Cord Digital LCD Display (LED based) Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Light Cure Unit With Cord Digital LCD Display (LED based) Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Light Cure Unit With Cord Digital LCD Display (LED based) China, Light Cure Unit With Cord Digital LCD Display (LED based) Kenya.

Light Cure Unit With Cord Digital LCD Display (LED based) Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Light Cure Unit With Cord Digital LCD Display (LED based) Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada

(Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point), Suriname (Paramaribo)

Technical Specification :

Light Cure Unit With Cord Digital LCD Display (LED based) Curing handpiece. Can be fixed on the chair. The special high performance diode mechanism can polymerize an average of 3mm thickness of material within 10 seconds. The smart three setting: Full Power, Soft-up and Impluse LED lamp 5 Watt blue lighting diode. LED lamp life time up to - 100, 000hrs. Wavelength - 450-490nm. Luminosity - 800-1200mw/cm

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [Medical Hospital Equipments Manufacturers](#),