

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Luxury Electric Patient Chair

**Product Code :**  
ENT0003



**Description :**

Luxury Electric Patient Chair The seat can be up & down, backrest can be forward and backward, arms moves together with leg pads, seat can be rotated 180 degree. Voltage : AC220V 50 HZ Power : 230 W (max) and low power consumption, energy saving. Under weight pressure : (minimum) 530mm – 680mm (maximum) Armrest and backrest connections moving, the backrest can be forwarded 90° into a line Pillow flex : 100mm (folding, detachable) Seat can be rotated 180°

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)