



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
SST (Gel Tubes) (MOQ: 100,000 Pcs Assorted Size)

**Product Code :**  
Blood0002



**Description :**

SST (Gel Tubes) (MOQ: 100,000 Pcs Assorted Size) Serum Separating Tubes(SSTs) contain the clot activator and a special gel at the base. The gel forms a stable barrier between the serum and the blood cells.After centrifugation, the serum sample can be used for special diagnostics.Mixing instructions: SSTs should be gently inverted 180° back and forth 5-6 times.Centrifugation instructions: 1300g-2000g for 10 minutes at 25° C room temperature



**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,