



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
SST (Gel Tubes) (MOQ: 100,000 Pcs Assorted Size)

Product Code :
Blood0002



Description :

SST (Gel Tubes) (MOQ: 100,000 Pcs Assorted Size) Serum Separating Tubes(SSTs) contain the clot activator and a special gel at the base. The gel forms a stable barrier between the serum and the blood cells. After centrifugation, the serum sample can be used for special diagnostics. Mixing instructions: SSTs should be gently inverted 180° back and forth 5-6 times. Centrifugation instructions: 1300g-2000g for 10 minutes at 25° C room temperature



Medical Hospital Equipments Manufacturers

Website: [SITE](#)