

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Sodium Citrate Tubes (MOQ: 100,000 Pcs Assorted Size)

Product Code :
Blood0004



Description :

Sodium Citrate Tubes (MOQ: 100,000 Pcs Assorted Size) Sodium Citrate Tubes blood collection tubes contain a measured amount of citrate solution and are used for coagulation determinations in citrated plasma. Since the liquid citrate dilutes the blood, it is essential that the tube is completely filled with blood sample so that the dilution is properly carried out. NET Sodium Citrate tubes ensure a constant mixing ratio of one part citrate to nine parts of blood

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [MedicalHospitalEquipmentsManufacturers,](http://MedicalHospitalEquipmentsManufacturers.com)