



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Multi - Sample Blood Collection Needles With Hub (MOQ:  
100,000 Pcs Assorted Size)

**Product Code :**  
Blood0007



**Description :**

Multi - Sample Blood Collection Needles With Hub (MOQ: 100,000 Pcs Assorted Size) The NET Multi -Sample Blood Collection Needles With Hub with optical venipuncture control can make an important contribution to simplifying blood collection routine for phlebotomists. The needle has a transparent plastic hub, which functions as a view window. If venipuncture is successful, the view window is coloured red due to the blood flow



**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,