

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Leur Adapter (MOQ: 100,000 Pcs Assorted Size)

**Product Code :**  
Blood0010



**Description :**

Leur Adapter (MOQ: 100,000 Pcs Assorted Size) The NET Leur adapter is a sterile ,single use device to be used with safety Blood collection sets or luer needles in routine blood collection .It is manufactured from stainless steel and polypropylene ans is fitted with a safety valve at one end.The perforated label not only serves to simplify identification ,but also acts as a seal if integrity

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,