

**Product Name :**  
Digital Portable Ultrasound with Led Monitor**Product Code :**  
NMCBS28003**Description :**

Digital Portable Ultrasound with Led Monitor Applications : Abdomen, Obsterics, Gynecology, Urology, Vascular, Small Parts, Pediatrics, Cardiology. Display Modes : B, B/B, 4B, B/M Gray Scale : 256 Levels Scanning Depth : 240mm Max (depend on the probe types) Probe Frequency : 2.5MHz – 11MHz (Depend on the probe types) Cine Loop : 256 frames Software Packages : General measurement, OB, Obstetrics, GYN, Urology, Abdomen, Cardiac B/M mode, Vascular, Small Part, Orthopedic, Pediatrics. Archives Management : Patient information management, Patient images record, Patient exam report. Documentation : Full screen documentation, Hospital name, Patient name, Date / time Annotation database available, Pre-saved phrases can be recalled Power Supply : AC 20V – 100V, 50Hz – 60Hz Standard Configuration : Main Unit + C3 – A Convex (3.5MHz) probe + 1 probe connector + 3 USB ports + VGA port + Video port Optional : 2 probe connectors, Built-in Battery, DICOM 3.0, Anti-water keyboard cover, Carry bag, Support PC Printer