

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Mammography System with Fixed Anode

**Product Code :**  
NMCBS28006



**Description :**

Mammography System with Fixed Anode Inverter frequency : 80kHz Input Power : 220VAC, 50/60Hz (can change to 110V) Power Rating : 3kW Focal : Small Focal Point 20-35kV/10-100mAs X-ray Tube : Fixed anode, XD7 Focal Spot : Dual Focus 0.5mm CIR : 18x24cm Weight : 248kg(1.58m3) Radiographic Stand : C-ARM : Vertical Movement, 580mm; Rotations Degree, +175° ~ 180°, SID : 650mm Anode : Fixed

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,