

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Mammography System with Rotating Anode

Product Code :
NMCBS28007



Description :

Mammography System with Rotating Anode Inverter frequency : 80kHz Input Power : 220VAC, 50/60Hz (can change to 110V) Power Rating : 6kW Focal : Large Focal Point 20-38kV/10-510mAs, Small Focal Point 20-35kV/10-100mAs X-ray Tube : LR01 Tube Focal Spot : Dual Focus 0.2 / 0.4mm Cassette Image Receptor : 18x24cm Weight : 300kg(1.68m³) Radiographic Stand : C-ARM : Vertical Movement : 580mm, Rotations Degree : +175° ~ 180°, SID : 650mm Anode : Rotating

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,