

Product Name :
Mobile X-Ray Machine**Product Code :**
NMCBS28009**Description :**

Mobile X-Ray Machine Use : NET mobile x ray equipments used in operation rooms, wards and emergency rooms etc. It serves the radiographic requirement of head, chest and limbs. Features : Flexible movement with Compact appearance, convenient in operation. With high voltage electrical tube with enough power. Three modes (15 mA, 30mA, 50mA), Different current can be selected as per the photography requirement. Anti electric shock, anti scattering and full wave rectifier. Stable and reliable electric performance. 120° Telescopic rotating angle

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,